


Notch rejection filters

# TopNotch™


Chroma Technology manufactures high-performance optical interference filters using sputtering and ion beam technologies. Covering a spectral range from 200nm- 3000nm, our wide-ranging capabilities enable us to offer a broad array of filter products offering superior durability and longevity, and spanning many industrial applications.




## Best-in-Class Performance

- Extra-wide, high-transmission range:  
350 - 1600nm
- Narrow notch:  
approx. 3% of notch CWL
- High blocking level:  
>OD6
- Wide blocking range:  
>6nm wavelength range

**Narrow band, square rejection**



**Wider range of high, OD6 blocking**



ZET488 TopNotch
  ZET532 TopNotch
  ZET594 TopNotch
  ZET642 TopNotch
  ZET785 TopNotch

